

Riduttori Per Provette - Nesting Sample cup




BSA636

Sample Nest cup Vol. 1ml riduttore per provette Ø 13 mm, in polistirolo.

Sample Nest Cup Vol.1 ml reduction cups, Ø 13 mm.

H 2,97 Ø 0,94 (max Ø 1,38) cm



BSA646

Sample Nest cup Vol. 2ml riduttore per provette Ø 16 mm, in polistirolo.

Sample Nest Cup Vol.2 ml reduction cups, Ø 16 mm.

H 2,97 Ø 1,24 (max 1,67) cm