

Provette Sorwall - Sorwall Test-Tubes



BSA059

Vol. 5ml, Ø 12x75 mm in polistirolo cilindrica.
Vol. 5ml, Ø 12x75 mm in polystyrene, cylindrical.

BSA059/P

Vol. 5ml, Ø 12x75 mm in polipropilene cilindrica.
Vol. 5 ml, Ø.12 X 75 mm in polypropylene, cylindrical.



BSA0592

Provetta Sorwall Ø12x75mm in PS blu.
Blue Sorwall test-tube in PS