

Provette Sterili - Sterilized Biotube

Provette per anaerobiosi con 2 punti di chiusura del tappo - Two-position closure cap Test-Tube




BSF2054

Provette graduate Vol. 5 ml, Ø 12 x 75 mm in polistirolo trasparente, sterili-
Test-tubes Vol. 5 ml graduated, Ø 12 x 75 mm in polystyrene sterile.



BSF2051

Provette graduate Vol. 14 ml, Ø 17 x 100 mm in polistirolo trasparente, sterili
Test-tubes Vol. 14 ml graduated, Ø 17 x 100 mm in polystyrene sterile.



Provette con tappo a vite sterili in polistirolo trasparente - Screw cap Test-Tube in polystyrene, sterile



BSM470

Provette tappo a vite vol. 10 ml, Ø 16 x 100 mm in polistirolo trasparente, sterili, fondo cilindrico.
Screw cap test-tube vol. 10 ml, Ø 16 x 100 mm in polystyrene sterile, cylindrical.

BSM470/CS

Idem BSM470 in confezione singola.
As BSM470 individually wrapped.



BSP250

Provette tappo a vite vol. 10 ml, Ø 16 x 100 mm in polistirolo trasparente, trattate, fondo conico. Disponibili anche con tappo rosso e bianco.

*Test-tubes screw cap vol. 10 ml, Ø 16 x 100 mm in polystyrene treated, conical
Also available with red and white stopper.*

Provette tappo a pressione sterili in polistirolo trasparente - Pressure Cap test-tubes sterile in polystyrene



BSM416

Provetta vol. 5 ml Ø 12 x 86 mm sterile, graduate 2,5 / 5 ml tappate, cilindriche.
Test-tubes vol. 5 ml Ø 12 x 86 mm sterile, graduated at 2,5 / 5 ml with stopper, cylindrical.



BSM400

Provetta vol. 10 ml Ø 16 x 100 mm sterile, graduate 2,5 / 5 / 10 ml tappate, cilindriche.
Test-tubes vol. 10 ml Ø 16 x 100 mm sterile, graduated at 2,5 / 5 / 10 ml with stopper, cylindrical.

BSM402

Idem BSM400 etichettate.
As BSM400 labelled.



BSM408

Provetta vol. 10 ml Ø 16 x 103 mm sterile graduate 1 / 2,5 / 5 / 10 ml tappate, coniche.
Test-tubes vol. 10 ml Ø 16 x 103 mm sterile graduated at 1 / 2,5 / 5 / 10 ml with stopper, conical.

BSM410

Idem BSM408 etichettate.
As BSM408 labelled.



BSM424

Provetta vol. 20 ml Ø 16 x 150 mm sterile, tappate, cilindriche.
Test-Tubes vol. 20 ml Ø 16 x 150 mm sterile, with stopper, cylindrical.

Provette Sterili con tappo a vite in polipropilene - *Sterile polypropylene test-tube, screw cap*

Provette sterili da 15 e 50 ml:

Serigrafia e banda di scrittura ben evidente e antigraffio, tappo a vite con chiusura ermetica.

Sterile centrifuge tubes 15 and 50 ml:

Permanent and non-scratch printed graduation, printed writing area, screw cap leak+ proof.



BSM482

Provette tappo a vite vol. 15 ml, \varnothing 17 x 120 mm in pp, graduate sterili a sterile fondo conico.

Test-tubes screw cap vol. 15 ml, \varnothing 17 x 120 mm in pp, graduated sterile, conical.

BSM490

Idem BSM482 in rack.

As BSM482 racked.



BSM474

Provette tappo a vite vol. 50 ml, \varnothing 30 x 115 mm in pp, graduate sterili a raggi beta fondo conico.

Test-tubes screw cap vol. 50 ml, \varnothing 30 x 115 mm in pp, graduated sterile, conical.

BSM475

Idem BSM474 in rack.

As BSM474 racked.



BSM477

Provette tappo a vite vol. 50 ml, \varnothing 30 x 115 mm in pp, graduate sterili fondo conico autportante.

Test-tubes screw cap vol. 50 ml, \varnothing 30 x 115 mm in pp, graduated sterile, self-standing.



BSM4701

Provette tappo a vite vol. 10 ml, \varnothing 16 x 100 mm in pp, sterili fondo cilindrico.

Test-tubes screw cap vol. 10 ml, \varnothing 16 x 100 mm in pp, sterile, cylindrical.

Provette tappo a vite in polistirolo trasparente, sterili - *Sterile screw cap Test-Tubes in polystyrene*



BSM483

Provette tappo a vite vol. 15 ml \varnothing 17 x 120 mm sterili, graduate, fondo conico in polistirolo.

Test-tubes screw cap vol. 15 ml \varnothing 17 x 120 mm sterile, graduated, conical in polystyrene.

BSM491

Idem BSM483 in rack.

As BSM483 in rack.

BSM485

Idem BSM483 in confezione singola.

As BSM483 individually wrapped.



BSM478

Provette tappo a vite vol. 16 ml \varnothing 16 x 120 mm sterili fondo cilindrico.

Test-tubes screw cap vol. 16 ml \varnothing 16 x 120 mm sterile cylindrical.

BSM478/CS

Idem BSM478 in confezione singola.

As BSM478 individually wrapped.



BSM486

Provette tappo a vite vol. 20 ml \varnothing 16 x 150 mm sterili fondo cilindrico.

Test-tubes screw cap vol. 20 ml \varnothing 16 x 150 mm sterile cylindrical.

BSM486/CS

Idem BSM486 in confezione singola.

As BSM486 individually wrapped.