

# Provette per Urina - Urine Biotube

Per M.T KIT e utilizzabili su strumenti **Superaution®**, **Supertron®**, **Atlas®** - Suitable to **Superaution®**, **Supertron®**, **Atlas®**



**BSP074**  
Provette per M.T. System® Vol. 10ml, 16x103 mm in polistirolo trasparente graduata a 1 / 2,5 / 5 / 10 ml.  
*M.T. System® Tube in polystyrene, Vol 10ml, 16x103 mm, graduated 1 / 2,5 / 5 / 10 ml*



**BSO026**  
Tappo Ø 16mm, colore bianco, in polietilene.  
*Stoppers Ø 16mm, white, in polyethylene*



**BSP100**  
Provette per test urina a fondo conico 12ml, in polistirolo.  
*Biotube Urine Test-Tube conical bottom 12 ml in polystyrene.*

**BSP107**  
Tappi per test urina fondo conico in polipropilene per BSP100.  
*Stoppers for test-urina conical bottom in polypropylene for BSP100.*



**BSP128**  
Provette per urina, tipo pediatrico, in polistirolo.  
*Test-Tube urine pediatric suitable to instruments, in polystyrene.*



**BSP075**  
Provette idem BSP074 etichettate.  
*Test-Tube as BSP074 labelled*



**BSP0741**  
Provette idem BSP074 con tappo a parte.  
*Test-Tube as BSP074, stopper aside*



**BSP111**  
Provette idem BSP074 tappate  
*Test-Tube as BSP074 capped*

**BSP112**  
Provette idem BSP074 tappate ed etichettate.  
*Test-Tube as BSP074 capped and labelled.*



**BSP113**  
Provette idem BSP074 in rack, tappate ed etichettate.  
*Test-Tube as BSP074 in rack capped and labelled.*