

Vasi Monouso - *Histology Containers*

Vaso monouso in polistirolo antiurto per pezzi chirurgici, con chiusura a pressione
Polystyrene disposable pot for surgical or anatomical specimens, with leak proof closure



Cod.	Dimension mm Ø	Volume ml	Serigraphy / Serigrafia
BSC2512	125x70	500	X
BSC2511	125x70	500	
BSC2522	125x120	1000	X
BSC2521	125x120	1000	
BSC2520	200x130	2500	X
BSC2524	200x130	2500	
BSC2527	200x215	5000	X
BSC2525	200x215	5000	