

Contenitori con Formalina - Containers with Formalin

Provette e contenitori in polipropilene. Tappo a vite, preriempita con formalina 4% neutra tamponata ph 7.2-7.4 pronto uso.
Prefilled tube / Containers, screw cap, with buffered 4% formalin ready to use.

REF.	TIPO	TYPE	DIM. CM	VOL. ML	VOL. RIEMPIMENTO ML / FILLING VOLUME ML
BSCQ01000	Provetta	Tube	H 13,50 Ø 1,60	10	5
BSCQ01011	Contenitore	Container	H 6,95 Ø 3,30	30	10
BSCQ01012	Contenitore	Container	H 6,95 Ø 3,30	30	20
BSCQ01021	Contenitore	Container	H 6,95 Ø 3,95	60	20
BSCQ01022	Contenitore	Container	H 6,95 Ø 3,95	60	40
BSCQ01023	Contenitore	Container	H 6,95 Ø 3,95	60	30
BSCQ01101	Contenitore	Container	H 10,25 Ø 6,05	150	90
BSCQ01111	Contenitore	Container	H 9,95 Ø 6,50	250	150



Contenitori in polipropilene. Tappo a pressione, preriempita con formalina 4% neutra tamponata ph 7.2-7.4 pronto uso.
Containers pressure cap, with buffered 4% formalin ready to use.

REF.	TIPO	TYPE	DIM. CM	VOL. ML	VOL. RIEMPIMENTO ML / FILLING VOLUME ML
BSCQ01061	Contenitore	Container	H 10,20 Ø 10,95	500	300
BSCQ01071	Contenitore	Container	H 13,20 Ø 13,50	1000	600
BSCQ01081	Contenitore	Container	H 14,50 Ø 19,50	2500	1500
BSCQ01091	Contenitore	Container	H 19,50 Ø 22,50	5000	3000